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Aniline blue staining experimental report

1. Experimental equipment and reagents

1.1 Experimental equipment

Name	Manufacturer	Model
CO2 incubator	SANYO	MCO-15AC
Clean Bench	Sujing Antai	SW-CJ-1FD
Slides	Wanwu	
Cover glass	Jiangsu Shitai Experimental Equipment Co.,Ltd	10127105P-G
ORTHO OPTICAL MICROSCOPE	NIKON, JAPAN	NIKON ECLIPSE E100
IMAGING SYSTEM	NIKON, JAPAN	NIKON DS-U3
Dehydrator	DIAPATH	Donatello
Embedding machine	Wuhan Junjie Electronics Co., Ltd	JB-P5
Pathology slicer	Shanghai Leica Instrument Co., Ltd	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd	JB-L5
Tissue spreader	Zhejiang Kehua Instrument Co., Ltd	KD-P
Oven	Tianjin Laibo Rui Instrument Equipment Co.,Ltd	GFL-230

1.2 Main experimental reagents

Reagent name	Manufacturer	Article number
Washing liquid	Wanwu	
Aniline blue staining solution	Wanwu	
Decolorizing solution	Serbicebio	
Neutral resin	Sinopharm Group Chemical Reagent Co. Ltd	10004160

2. Experimental

For plants

2.1 The decolorized leaves were taken, embedded and sliced, then dewaxed to water and washed with pure water for three times.

2.2 Wash the slides with detergent for 3 times, dry the detergent and put into 0.01% aniline blue dye, dip away from light at room temperature for 2 hours.

2.3 Seal: The slides were washed with detergent for 3 times to absorb dry water. Drop anti-fluorescence quenching sealing film liquid sealing sheet, cover glass slide.

2.4 Take photos: observe and take photos. Observe the deposition of callose with fluorescence microscope (DAPI filter). The fluorescence is easy to quench. Take photos as soon as possible and avoid light. The slices can be stored at 4 °C for a short time.

For blades

2.1. Decolorization: fresh plant tissues (young leaves) are decolorized with decolorizing solution.

2.2 Dyeing: the tissue was washed with washing solution for 3~5 times, the excess liquid was dried by filter paper, and bleached in 0.01% aniline blue dyeing solution at room temperature without light for 2 hours.

2.3 Spreading and sealing: the tissue is rinsed with washing solution for 3~5 times, the filter paper absorbs the excess liquid and lays it flat on the carrier, dripping with anti-fluorescent sealing agent, cover glass and gently flatten the tissue without air bubbles.

2.4 Take photos: observe and take photos. Observe the deposition of callose with fluorescence microscope (DAPI filter). The fluorescence is easy to quench. Take photos as soon as possible and avoid light. The slices can be stored at 4 °C for a short time.

3. The results were as follows

3.1 The callose showed bright blue fluorescence and was widely distributed around the sieve tub.