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## Luxol Fast Blue(LFB) staining experimental report

### 1. Lab equipment and reagents

#### A. Lab equipment

Items	Manufacturer	Model
Dehydrator	DIAPATH	Donatello
embedding machine	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathology microtome	Shanghai Leica Instruments Co., Ltd.	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Water Bath-Slide Drier	Zhejiang Jinhua Kedi Instrumental Equipment CO.,LTD	KD-P
Laboratory oven	Tianjin Labotery Instrument Equipment Co., Ltd.	GFL-230
Microscope slide	Wanwu	
Upright optical microscope	Nikon Japan	Nikon Eclipse E100
Imaging system	Nikon Japan	NIKON DS-U3

#### B. Chemical Reagents

Items	Manufacturer	Model
Absolute alcohol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
Luxol fast blue staining kit	Wanwu	G1030
Hematoxylin-Eosin solution	Wanwu	G1002
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

### 2. Experimental steps

(1) Paraffin section deparaffinization and rehydration: put the slides into xylene I 20 minutes-xylene II 20 minutes-absolute ethanol I 5 min-absolute ethanol II 5 min-75% alcohol for 5 min, tap water washing.

(2) Luxol fast bluestaining: placed luxol fast blue staining A is in a 60°C oven and preheated for 30 minutes. Put he slides intoluxol fast blue staining A and cover with membrane for 1 hour, then take out the slides quickly wash with tap water.

(3) Background differentiation: put the slides into luxol fast blue staining B for 2s and directly immerse in luxol fast blue staining C for differentiation 15s, washing to terminate differentiation, microscope examination, repeated differentiation, washing and microscopic examination, until the myelin sheath is blue, the background is almost colorless;

(4) Eosin counterstaining: put the slides in the oven at 65 degrees to dry it (about 30mins until the slides are dry) and take out and cool it, put into 95% ethanol and Eosin counterstaining.

(5) Dehydration and sealing: put the slides into absolute ethanol I 5min-anhydrous ethanol II 5min-anhydrous ethanol III 5min-xylene I 5min-xylene II 5min -transparent, neutral balsam sealing.

(6) Microscope examination, images collection and analysis.

### 3. Results

Nerve myelin sheath is blue, and other components are almost colorless.

### 4. Note

1. After stained, do not wash them in water for a long time due to the color will be fade.
2. Pay attention to the degree of alternating differentiation, if over-differentiation, the color of the myelin sheath will be fade, insufficient differentiation, the boundaries are unclear.
3. If over-differentiation, you can re-stain the slides in the staining solution.