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β -galactosidase Staining Experiment Report

I. Experimental equipment and reagents

1. Experimental equipment

| Equipment name | Manufacturer | Model No. |
|-----------------------------|-------------------------------------|--------------------|
| Freezing microtome | Thermo | CRYOSTAR NX50 |
| Slicing knife | Shanghai Leica Instrument Co., Ltd. | LEICA 819 |
| Microscope slides | Wanwu | |
| Upright electron microscope | Japan NIKON | NIKON ECLIPSE E100 |
| Imaging system | Japan NIKON | NIKON DS-U3 |

2. Main experimental reagents

| Reagent name | Manufacturer | Item No. |
|--|--------------|----------|
| OCT Embedding medium | Sakura | 4583 |
| β -galactosidase staining solution A | Wanwu | G1073 |
| β -galactosidase staining solution B | Wanwu | G1073 |
| β -galactosidase staining solution C | Wanwu | G1073 |
| β -galactosidase staining solution D | Wanwu | G1073 |
| β -galactosidase staining solution E | Wanwu | G1073 |
| β -galactosidase staining solution F | Wanwu | G1073 |
| Glycerin Jelly | Wanwu | G1402 |

II. Staining steps

1. Preparation of working solution: 940ul of staining solution A, 10ul of staining solution B, and 50ul of β -galactosidase staining solution D, add them in sequence and mix well, then the working

solution is made.

2. Fixation: Take a fresh frozen section that has not been fixed and let it dry naturally, fix it with β -galactosidase staining solution E for 15 min, pour out the fixation solution, and wash it three times with distilled water.

3. β -galactosidase staining: Spin off the water on the section, drop the working solution onto the tissue, and incubate at 37°C for 16 to 18 hours.

4. Pour out the working solution and wash adequately with distilled water for 3 times.

5. After incubation with β -galactosidase stain F for 3 min, wash with water.

6. Dehydration and sealing: Put the sections in sequence into Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - Anhydrous ethanol III for 5min - Xylene I for 5min - Xylene II for 5min for transparency, seal with neutral gum lastly.

III. Interpretation of results:

Senescent cells are blue.

IV. Precautions:

1. Specimens need to be fresh frozen sections or cells, and can't be fixed.
2. The staining solutions need to be stored at -20°C, and the shelf life is one year.
3. The working solution should be prepared freshly and use immediately, and cannot be stored.