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1.2 Major reagents

Name	Producer	Code
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
Masson dye solution set	Wanwu	G1006
Neutral gum	SCRC	10004160

2 Procedure

2.1 Dewaxing as followed:

- Xylene I for 20 min;
- Xylene II for 20 min;
- 100% ethanol I for 5 min;
- 100% ethanol II for 5 min;
- 75% ethanol for 5 min;
- Rinsing with tap water ;

- 2.2 Soak the slices in Masson A overnight, rinse with tap water.
- 2.3 Masson B, Masson C were prepared into Masson solution according to the ratio of 1:1. Then stain with Masson solution for 1 min, rinse with tap water. And differentiation it with 1% hydrochloric acid alcohol, rinse with tap water.
- 2.4 Soak the slices in Masson D for 6 min, rinse with tap water; Masson E for 1 min; Masson F for 2-30s.
- 2.5 Rinse the slices with 1% glacial acetic acid and then dehydration with two cup of anhydrous ethanol.
- 2.6 100% ethanol for 5 min; Xylene for 5 min; Finally seal with neutral gum.
- 2.7 Observe with microscope inspection, image acquisition and analysis.

3 Results

Color	Result
Blue	Collagen fibers
Red	Muscle fibers; Cellulose; Erythrocyte

4 Precautions

- 4.1 After treat the slices with Masson A , raine with the water should be quickly. It is sufficient to just wash the slices clean.
- 4.2 Replace the mixed dye solution according to the number of slices. And it can also be formulated when you need.
- 4.3 The staining degree of Masson D solution should be well controlled. The collagen part cannot be red, which will affect the subsequent coloration of Masson F solution;
- 4.4 During the differentiation of glacial acetic acid, if it is over-differentiated, the collagen blue will be too light, if it is insufficiently differentiated, it will be easily superimposed with red to become purple-blue .