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## **VG staining experimental report STAINING- EXPERIMENT REPORT**

### **1. Experimental equipment and reagents**

#### **1.1 Experimental equipment**

<b><u>Name</u></b>	<b><u>Manufacturer</u></b>	<b><u>Model</u></b>
<u>Dehydrator</u>	<u>DIAPATH</u>	<u>Donatello</u>
<u>Embedding machine</u>	<u>Wuhan Junjie Electronics Co., Ltd</u>	<u>JB-P5</u>
<u>Pathology slicer</u>	<u>Shanghai Leica Instrument Co., Ltd</u>	<u>RM2016</u>
<u>Frozen platform</u>	<u>Wuhan Junjie Electronics Co., Ltd</u>	<u>JB-L5</u>
<u>Tissue spreader</u>	<u>Zhejiang Kehua Instrument Co., Ltd</u>	<u>KD-P</u>
<u>Oven</u>	<u>Tianjin Laibo Rui Instrument Equipment Co., Ltd</u>	<u>GFL-230</u>
<u>Slides</u>	<u>Wanwu</u>	
<u>Orthostatic microscope</u>	<u>NIKON, JAPAN</u>	<u>NIKON ECLIPSE E100</u>
<u>Image system</u>	<u>NIKON, JAPAN</u>	<u>NIKON DS-U3</u>

### **1. Experimental equipment and reagents**

experiment-equipment

<b><u>name</u></b>	<b><u>factory</u></b>	<b><u>model</u></b>
<u>Dehydrator</u>	<u>DIAPATH</u>	<u>Donatello</u>
<u>Embedding machine</u>	<u>Wuhan Junjie Electronics Co., Ltd.</u>	<u>JB-P5</u>
<u>Pathology slicer</u>	<u>Shanghai Leica Instruments Co., Ltd.</u>	<u>RM2016</u>
<u>Frozen platform</u>	<u>Wuhan Junjie Electronics Co., Ltd.</u>	<u>JB-L5</u>
<u>Organizer</u>	<u>Jinhua City, Zhejiang Province Cody</u> <u>Instrument-Equipment-Co.,Ltd.</u>	<u>KD-P</u>
<u>oven</u>	<u>Tianjin=Lai=Bo=Rui=Instrument=Equipment=GFL-230</u> <u>Co.,Ltd.</u>	
<u>Glass slide</u>	<u>Wanwu</u>	
<u>UPRIGHT OPTICAL</u> <u>MICROSCOPE</u>	<u>NIKON JAPAN</u>	<u>NIKON ECLIPSE E100</u>

IMAGING SYSTEM

NIKON JAPAN

NIKON DS-U3

## 1.2 Main experimental reagents

### Main experimental reagents

<u>Reagent name</u>	<u>Reagent name</u>	<u>Manufacturer</u>	<u>factory</u>	<u>Article number</u>	<u>Article number</u>
Absolute ethanol		Sinopharm Group Chemical Reagent Co., Ltd.		100092683	
Xylene		Sinopharm Group Chemical Reagent Co., Ltd.		10023418	
VG <u>staining</u>		Wanwu		G1046	
Neutral gum		Sinopharm Group Chemical Reagent Co., Ltd.		10004160	

## 2. Experimental steps

2.1 Paraffin slides dewaxed as follow: Two changes of pure xylene for 20 min; Two changes of pure ethanol for 5min; 75% ethanol for 5min; Keep slides in tap water.

2.2 VG dyeing: Add 9ml of Van Gieson (VG) dye solution B to 1ml of Van Gieson (VG) dye solution A (can be used repeatedly after preparation) and mix into VG dye solution (prepared in proportion to use as much as possible), dye for 1 minute, fast Washed with water, dehydrated quickly in anhydrous ethanol triplex.

2.3 ~~III. The slides were immersed in xylene to transparent and then coverslip with neutral resin. Transparent sealing: two cylinders of clean xylene are transparent for 20s and 5min each (xylene is not shared with other xylenes), neutral gum wet seal.~~

~~IV 2.4. Observed under microscope and took images.~~

~~Microscope inspection, image acquisition and analysis.~~

## 3. The results were as follows

Collagen fibers ~~was~~ red, and other tissue components ~~was~~ yellow.

## 4 Four, Precautions

1. It should be washed quickly to prevent the red from fading, and the absolute ethanol should also be taken off quickly to prevent the yellow from darkening.