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1.2 Major reagents

Name	Producer	Code
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
Congo Red dye solution set	Wanwu	G1003
Differentiation solution	Wanwu	G1005-3
Scott Tap Bluing	Wanwu	G1005-4
Neutral gum	SCRC	10004160

2 Procedure

2.1 Dewaxing as followed:

- Xylene I for 20 min;
- Xylene II for 20 min;
- 100% ethanol I for 5 min;
- 100% ethanol II for 5 min;
- 75% ethanol for 5 min;
- Rinsing with tap water ;

2.2 Soak the slices in Congo Red A overnight, rinse with tap water for 2 min.

- 2.3 Treat the section with Congo Red B for 1s until the positive plaques were obvious, and the background was basically colorless, then rinse with tap water.
- 2.4 Stain sections with Congo Red C for 1 min, rinse with tap water. Then treat it with Differentiation solution and Scott Tap Bluing, each step required washing with water.
- 2.5 Dehydrate as followed:
 - 100% ethanol I for 5 min;
 - 100% ethanol II for 5 min;
 - 100% ethanol III for 5 min;
 - Xylene I for 5 min;
 - Xylene II for 5 min;
 - Finally seal with neutral gum.
- 2.6 Observe with microscope inspection, image acquisition and analysis.

3 Results

Color	Result
Blue	Nucleus
Red (under white light)	Amyloid

4 Precautions

- 4.1 Please pay attention to distinguishing thyroid glial, elastic fiber and collagen fiber, which may be dyed red together;
- 4.2 It is necessary to master proper differentiation, insufficient differentiation, collagen fibers will also become red, over-differentiated, amyloid protein will also be decolorized.
- 4.3 If there is no positive or positive weakly under the white light, you can observe the red light under fluorescence.