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Safranin fast green staining for bone experimental report

1. Experimental equipment and reagents

1.1 Experimental equipment

| Name | Manufacturer | Model |
|------------------------------|--|--------------------|
| Dehydrator | Wuhan Junjie Electronics Co., Ltd | JJ-12J |
| Embedding machine | Wuhan Junjie Electronics Co., Ltd | JB-P5 |
| Pathological section machine | Shanghai Leica Instrument Co., Ltd | RM2016 |
| Frozen platform | Wuhan Junjie Electronics Co., Ltd | JB-L5 |
| KD-P Water Bath | Kedee | KD-P |
| Oven | Tianjin Lai Bo Rui Instrument Equipment Co., Ltd | GFL-230 |
| Slides | Wanwu | |
| Orthostatic microscope | NIKON, JAPAN | NIKON ECLIPSE E100 |
| Image system | NIKON, JAPAN | NIKON DS-U3 |

1.2 Main experimental reagents

| Reagent name | Manufacturer | Article number |
|--|--|----------------|
| Absolute ethanol | Sinopharm Group Chemical Reagent Co. LTD | 100092683 |
| Xylene | Sinopharm Group Chemical Reagent Co. LTD | 10023418 |
| Saffron Fast Green dye solution (Bone Tissue) | Wanwukit | G1053 |
| Differentiation solution | Wanwu | G1040 |
| Neutral resin | Sinopharm Group Chemical Reagent Co. LTD | 10004160 |

2. Experimental steps

2.1 Paraffin slides dewaxed as follow: Two changes of pure xylene for 20 min. Two changes of pure ethanol for 5min. 75% ethanol for 5min. Keep slides in tap water.

2.2 Fast green staining: The slides stained in fast green dye solution for 1-5 minutes, washed away the excess dye solution until the cartilage was colorless, and soaked in 1% hydrochloric acid and alcohol for 10s. And then washing them with tap water.

2.3 Saffron staining: The slides were stained in saffron dye solution for 1-5s, and then put into four cylinders of absolute ethanol, where three for rapid dehydration for 5s, 2s, 10s respectively, and kept in the fourth cylinder.

2.4 The slides were immersed in xylene to transparent for 5min, sealing with neutral resin.

2.5 Observed under microscope, and took images for analysis.

3. The results were as follows

Cartilage was red or orange-red, and bone formation was green. And some connective tissues were red.

4. Precautions

4.1 During the dyeing process, the saffron staining cannot be overtime, because it is easy to hybridize with fast green to become purple-blue.

4.2 If the result of the staining is not well, it can be washed with hydrochloric acid solution and alcohol until the color fades, and then re-stained.