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## Copper red acid staining experimental report

### 1. Lab equipment and reagents

#### A. Lab equipment

Items	Manufacturer	Model
Dehydrator	DIAPATH	Donatello
embedding machine	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathology microtome	Shanghai Leica Instruments Co., Ltd.	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Water Bath-Slide Drier	Zhejiang Jinhua Kedi Instrumental Equipment CO.,LTD	KD-P
Laboratory oven	Tianjin Labotery Instrument Equipment Co., Ltd.	GFL-230
Microscope slide	Wanwu	
Upright optical microscope	Nikon Japan	Nikon Eclipse E100
Imaging system	Nikon Japan	NIKON DS-U3

#### B. Chemical Reagents

Items	Manufacturer	Model
Absolute alcohol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
Dithio oxamide	Sigma	379387
Acetic acid	Sinopharm Chemical Reagent Co., Ltd.	10018818
Nuclear fast red	Wanwu	G1035
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

### 2. Reagent preparation

0.1% Dithio oxamide ethanol solution: Dithio oxamide 10mg + absolute ethanol 10ml

10% acetic acid: acetic acid 10g + ultra-pure water 100ml

### 3. Experimental steps

- 
- (1) Paraffin section deparaffinization and rehydration: put the slides into xylene I 20minutes-xylene II 20 minutes-absolute ethanol I 5 min-absolute ethanol II 5 min-75% alcohol for 5 min, then tap water washing and distilled water washing for 3 times.
  - (2) Copper red acid staining: 0.1% Dithio oxamide ethanol solution: 10% sodium acetate mix at a ratio of 1:20, and put the slides into the mixed solution and covered in a 37°C water bath for about 3 days. Put the slides into 70% ethanol twice for 10 minutes each time, and absolute ethanol washing for twice, 3 hours each time;
  - (3) Nuclear staining: distilled water washing for 3 times, wash out ethanol, nuclear fast red counterstaining for 3min, then tap water washing;
  - (4) Dehydration and sealing: dehydrate with absolute ethanol in 3 cylinders, transparent with n-butanol for 5min, and xylene for 5min, neutral balsam sealing.
  - (5) Microscope examination, image collection and analysis.

#### 4. Results

The copper salt deposition are dark green and black, and the nucleus is red.

#### 5. Note

1. 0.1% Dithio oxamide ethanol solution should be stored at room temperature( lightproof);
2. The copper red acid mixed solution should be ready for use and discard after usage.