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#### 1.2 Major reagents

<b>Name</b>	<b>Producer</b>	<b>Code</b>
Xylene	SCRC	10023418
Toluidine Blue	Wanwu	G1032
Glacial acetic acid	Wanwu	G10000218
Neutral gum	SCRC	10004160

### **2 Procedure**

#### 2.1 Dewaxing as followed:

- Xylene I for 20 min;
- Xylene II for 20 min;
- 100% ethanol I for 5 min;
- 100% ethanol II for 5 min;
- 75% ethanol for 5 min;
- Rinsing with tap water ;

2.2 About plant: Treat the plant tissue slices with Toluidine Blue for 2-5 minutes, rinse with tap water. Then microscopically inspected it, and according to the color of the tissues, whether they are needed differentiation. Wash it with tap water, then dry it in the oven.

2.3 About animal: Treat the animal tissue slices with Toluidine Blue for 2-5 minutes, rinse with tap water. Treat it with 1% Glacial acetic acid. And you can use tap water to stop the

differentiation. Control the degree of differentiation under a microscope. Wash it with tap water, then dry it in the oven.

2.4 Xylene for 10 min, and seal with neutral gum. Please do not use xylene used in other experiments.

2.5 Observe with microscope inspection, image acquisition and analysis.

### 3 Results

Project	Color	Result
Plant	Blue-green	Lignified cell wall
	Purple-blue	Cellulose cell wall
Brain	Navy blue	Nissl
	Light blue	Background
Animal Bone	Purple blue	Cartilage
	Light blue	Background
Intestinal	Fuchsia	Mast cells
	Light blue	Nucleus
	Light blue	Background

### 4 Precautions

4.1 After dyeing, the slices cannot be left in water for a long time, otherwise they will fade.

4.2 The tissue must be completely dried before sealing, to avoid produces water vapor.