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Bisulfite sequencing PCR (BSP) Lab Report

1 Laboratory equipment and reagents

1.1 Laboratory equipment

Equipment	Manufacturers	Model
Centrifuge	DragonLab	D3024R
PCR instrument	EASTWIN	ETC811
Constant temperature oscillation shaker	Changzhou Aohua Instrument Co., Ltd	SHZ-82A
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Electrophoresis	Wanwu	FW-600
Gel imaging system	Tanon Science & Technology Co., Ltd.	Tanon-1600R
Standard Reagent Type Lab Ultra pure Water Purifier	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Dry thermostat	Hangzhou Allsheng Instruments Co.,Ltd	K20
Ice maker	Changshou City shircore	IMS-20
Ultramicro spectrophotometer	Electrical Appliance Co. Ltd.	NanoDrop2000
Mini-centrifuge	Thermo	MC-700
Wanwu		
Vortex oscillator	Wanwu	MV-100
High-speed-microtherm		
Wanwu		
Homogenizer		KZ-III-F

1.2 Reagents

Reagents	Manufacturers	Order
TIANamp Genomic DNA Kit	Tiagen Biotech(Beijing) CO.,LTD	DP304
Methylation-Gold Kit	ZYMO	D5005S(10)
TIANGel Midi Purification Kit	Tiagen Biotech(Beijing) CO.,	DP209
TaKaRa Taq Hot Start DNA Polymerase	LTDTAKARA Bioiotechnology (Dalian)	DR007A
T-Vector PMD19	Co.,LTD	K1901AA
2×Fast Pfus PCR Master Mix	TAKARA	G3305
Agarose	Wanwu	G5056

Anhydrous ethanol	Sinopharm Group Chemical Reagent Co., Ltd.	10009218
Isopropyl alcohol	Sinopharm Group Chemical Reagent Co., Ltd.	80109218
MarkerI DNA ladder	Wanwu	WGM1100
Centrifuge tube(1.5 mL)	Wanwu	EP-150-M
TIP	Wanwu	
Gold view	Wanwu	G1066
Proteinase K	Wanwu	G1205
Red Blood Cell Lysis Buffer	Wanwu	G2015
50*TAE	Wanwu	G3001
10×TE (pH 8.0)	Wanwu	G3003
6× DNA Loading Buffer	Wanwu	G3011
Ethidium bromide	Wanwu	G3010
Luria-Bertani Liquid Medium (Powder)	Wanwu	G3102

Tube and TIP are moist heat sterilization 40min, dry。

2 Experimental steps

2.1 Genomic DNA Extraction Experimental steps

2.1.1 Samples preparation

2.1.1.1 For blood, please use 200 μ l fresh, frozen or anticoagulant adding blood. If less than 200 μ l, please make up with buffer GA to 200 μ l.

2.1.1.2 If the sample is blood from poultry, birds, amphibians, of which red blood cells have nucleolus, the amount should be reduced to 5-20 μ l and adjust the volume to 200 μ l with buffer GA.

2.1.1.3 The adherent cells should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ($\sim 11,200 \times g$), then discard the flow-through and re-suspend cell pellet in 200 μ l buffer GA.

2.1.1.4 Animal tissue (spleen < 10 mg) should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ($\sim 11,200 \times g$), then discard the flow-through and re-suspend cell pellet in 200 μ l buffer GA.

2.1.2 Add 20 μ l Proteinase K, mix thoroughly by vortex. If the sample is tissue: incubate at 56°C until the tissue is completely lysed.

2.1.3 Add 200 μ l Buffer GB to the sample, mix thoroughly by vortex, and incubate at 70°C for 10 min to yield a homogeneous solution. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.1.4 Add 200 μ l ethanol (96-100%) to the sample, and mix thoroughly by vortex for 15 s. A white precipitate may form on addition of ethanol. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.1.5 Pipet the mixture from step 4 into the Spin Column CB3 (in a 2 ml collection tube) and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 s. Discard flow-through and place the spin column into the collection tube.

2.1.6 Add 500 μ l Buffer GD (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 s, then discard the flow-through and place the spin column into the collection tube.

2.1.7 Add 600 μ l Buffer PW (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 s. Discard the flow-through and place the spin column into the collection tube.

2.1.8 Repeat Step 7.

2.1.9 Centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 2 min to dry the membrane completely. Note: The residual ethanol of buffer PW may have some affection in downstream application

2.1.10 Place the Spin Column CB3 in a new clean 1.5 ml microcentrifuge tube, and pipet 50-200 μ l Buffer TE directly to the center of the membrane. Incubate at room temperature (15-25°C) for 2-5 min, and then centrifuge for 2 min at 12,000 rpm ($\sim 13,400 \times g$).

2.2 Bisulfite treatment of DNA

2.2.1 Add 130 μ l of the CT Conversion Reagent to 20 μ l of your DNA sample in a PCR tube. If the volume of the DNA sample is less than 20 μ l, make up the difference with water. Mix the sample by flicking the tube or pipetting the sample up and down, then centrifuge the liquid to the bottom of the tube.

2.2.2 Place the sample tube in a thermal cycler and perform the following steps*:

1. 98°C for 10 minutes
2. 64°C for 2.5 hours
3. 4°C storage up to 20 hours.

2.2.3 Add 600 μ l of M-Binding Buffer to a Zymo-Spin™ IC Column and place the column into a provided Collection Tube.

2.2.4 Load the sample (from Step 2) into the Zymo-Spin™ IC Column containing the M-Binding Buffer. Close the cap and mix by inverting the column several times.

2.2.5 Centrifuge at full speed ($>10,000 \times g$) for 30 seconds. Discard the flow-through.

2.2.6 Add 100 μ l of M-Wash Buffer to the column. Centrifuge at full speed for 30 seconds.

2.2.7 Add 200 μ l of M-Desulphonation Buffer to the column and let stand at room temperature (20-30°C) for 15-20 minutes. After the incubation, centrifuge at full speed for 30 seconds.

2.2.8 Add 200 μ l of M-Wash Buffer to the column. Centrifuge at full speed for 30 seconds. Add another 200 μ l of M-Wash Buffer and centrifuge for an additional 30 seconds.

2.2.9 Place the column into a 1.5 ml microcentrifuge tube. Add 10 μ l of M-Elution Buffer directly to the column matrix. Centrifuge for 30 seconds at full speed to elute the DNA.

2.3 PCR amplification

2.3.1 0.2ml PCR tube, 20ul pcr amplification system.

10×PCR Buffer	5.0 μ l
2.5mM dNTPs	5.0 μ l
10 μ M primers(F+R)	2.0 μ l+ 2.0 μ l
DNA	1.5 μ l
5U/ μ l Taq HS	0.2 μ l
ddH ₂ O	34.3 μ l

2.3.2 PCR amplification

Pre-denaturation	95°C, 5min	} 40×cycle
Denatured	95°C, 30s	
Annealing	55°C, 30s	
Extended	72°C, 30s	
End extension	72°C, 5min	

2.4 Electrophoresis: 1×TAE, 2.0% agarose, 4V/cm

2.5 Purify DNA from agarose gel

2.5.1 Column equilibration: add 500 μ l Buffer BL to the Spin Column CA2 (put Spin Column CA2 into a collection tube). Centrifuge for 1 min at 12,000 rpm ($\sim 13,400 \times g$) in a table-top microcentrifuge. Discard the flow-through, and put Spin Column CA2 back into the collection tube.

2.5.2 Cut the DNA fragment from agarose gel with a clean, sharp scalpel. Weigh the gel slice in a clean tube.

2.5.3 Add equivalent volume of Buffer PN to the gel (If the gel is 0.1 g, it is defaulted to be 100 μ l, then add 100 μ l Buffer PN). Incubate at 50°C by inverting up and down the tube until the agarose gel dissolves completely. If the agarose gel does not dissolve completely, incubate for longer period or add additional Buffer PN until all the agarose gel dissolved completely.

2.5.4 When the gel dissolved completely and the solution temperature turns to room temperature (15-25°C), transfer the mixture to the Spin Column CA2 (put Spin Column CA2 into a collection tube). Let the column stand for 2 min at room temperature (15-25°C), then centrifuge for 30-60 s at 12,000 rpm ($\sim 13,400 \times g$) in a table-top microcentrifuge. Discard the flow-through; place the Spin Column CA2 back into the collection tube again.

2.5.5 Wash the Spin Column CA2 with 600 μ l Buffer PW (ensure that ethanol (96-100%) has been added) and centrifuge for 30-60 s at 12,000 rpm ($\sim 13,400 \times g$). Discard the flow-through and place the Spin Column CA2 back into the collection tube.

2.5.6 Repeat Step 5.

2.5.7 Place the Spin Column CA2 back to the collection tube and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 2 min to remove residual wash buffer. Discard the flow-through, and place column with the cap open for several minutes to air dry the membrane.

2.5.8 Transfer the Spin Column CA2 to a clean 1.5 ml microcentrifuge tube. Add appropriate volume of Buffer EB to the center of the membrane, incubate at room temperature (15-25°C) for 2 min, then centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 2 min.

2.6 Recommended DNA insert volume

2 $\times$ solution	6.0 μ l
T-Vector PMD19	1.0 μ l
purify DNA	5.0 μ l

16°C connection overnight.

2.7 Transformation experiment using DH5 α

2.7.1 Thaw competent cells on wet ice for 20 minutes.

2.7.2 For DNA from ligation reactions. Add 1 μ l of the dilution to the cells, moving the pipette through the cells while dispensing. Gently tap tubes to mix.

2.7.3 Incubate cells on ice for 20 minutes.

2.7.4 Heat-shock cells 90 seconds in a 42°C water bath; do not shake.

2.7.5 Place on ice for 2 minutes.

2.7.6 Add 0.8 ml of room temperature LB Medium.

2.7.7 Shake at 220 rpm (37°C) for 45 minutes.

2.8 PCR amplification identification

The positive colonies on the plate were inoculated into the liquid medium and shaken at 37 ° C.
pcr amplification for Positive clones.

2.8.1 0.2ml PCR tube, 20ul pcr amplification system.

2×Fast Pfus PCR Master Mix	25ul
10μM primers(F+R)	1.5μl+1.5μl
Bacteria liquid	2.0μl
ddH2O	20.0μl

2.8.2 PCR amplification

Pre-denaturation	98°C, 2min	} 30×cycle
Denatured	98°C, 20s	
Annealing	55°C, 20s	
Extended	72°C, 10s	
End extension	72°C, 5min	

2.9 DNA sequencing

DNA sequencing by Wuhan GeneCreate Biological Engineering Co., Ltd.