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2.4 电泳：1 $\times$ TAE，2.0% agarose，4V/cm

Methylation-specific PCR(MSP) Lab Report

1 Laboratory equipment and reagents

1.1 Laboratory equipment

Equipment	Manufacturers	Model
Centrifuge	DragonLab	D3024R
PCR instrument	Beijing Dongsheng Innovation Biotechnology Co., Ltd.	ETC811
Constant temperature oscillation shaker	Changzhou Aohua Instrument Co., Ltd	SHZ-82A
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Electrophoresis	Wanwu	FW-600
Gel imaging system	Tanon Science & Technology Co., Ltd.	Tanon-1600R
Standard reagent type pure water meter	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Dry thermostat	Hangzhou Allsheng Instruments Co.,Ltd	K20
Ice maker	Changshou City shircore Electrical Appliance Co. Ltd.	IMS-20
Ultramicro spectrophotometer	Thermo	NanoDrop2000
Mini-centrifuge	Wanwu	MC-700
Vortex oscillator	Wanwu	MV-100
High-speed-microtherm Homogenizer	Wanwu	KZ-III-F

1.2 Reagents

Reagents	Manufacturers	Order
TIANamp Genomic DNA Kit	Tiagen Biotech(Beijing) CO.,LTD	DP304
Methylation-Gold Kit	ZYMO	D5005S(10)
TaKaRa Taq Hot Start DNA Polymerase	TAKARA Bioiotechnology(Dalian) Co.,LTD	DR007A
Agarose	Wanwu	G5056
Anhydrous ethanol	Sinopharm Group Chemical Reagent Co., Ltd.	10009218
Isopropyl alcohol	Sinopharm Group Chemical Reagent Co., Ltd.	80109218
MarkerI DNA ladder	Wanwu.	WGM1100
Centrifuge tube(1.5 mL)	Wanwu	EP-150-M

TIP	Wanwu	
Gold view	Wanwu	G1066
Proteinase K	Wanwu	G1205
Red Blood Cell Lysis Buffer	Wanwu	G2015
50*TAE	Wanwu	G3001
10×TE (pH 8.0)	Wanwu	G3003
6× DNA Loading Buffer	Wanwu	G3011
Ethidium bromide	Wanwu	G3010
Luria-Bertani Liquid Medium (Powder)	Wanwu	G3102

2 Experimental steps

2.1 Genomic DNA Extraction Experimental steps

2.1.1 Samples preparation

2.1.1.1 For blood, please use 200 µl fresh, frozen or anticoagulant adding blood. If less than 200 µl, please make up with buffer GA to 200 µl.

2.1.1.2 If the sample is blood from poultry, birds, amphibians, of which red blood cells have nucleolus, the amount should be reduced to 5-20 µl and adjust the volume to 200 µl with buffer GA.

2.1.1.3 The adherent cells should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ($\sim 11,200 \times g$), then discard the flow-through and re-suspend cell pellet in 200 µl buffer GA.

2.1.1.4 Animal tissue (spleen < 10mg) should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ($\sim 11,200 \times g$), then discard the flow-through and re-suspend cell pellet in 200 µl buffer GA.

2.1.2 Add 20 µl Proteinase K, mix thoroughly by vortex. If the sample is tissue: incubate at 56°C until the tissue is completely lysed.

2.1.3 Add 200 µl Buffer GB to the sample, mix thoroughly by vortex, and incubate at 70°C for 10 min to yield a homogeneous solution. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.1.4 Add 200 µl ethanol (96-100%) to the sample, and mix thoroughly by vortex for 15 s. A white precipitate may form on addition of ethanol. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.1.5 Pipet the mixture from step 4 into the Spin Column CB3 (in a 2 ml collection tube) and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 s. Discard flow-through and place the spin column into the collection tube.

2.1.6 Add 500 µl Buffer GD (Ensure ethanol (96-100%) has been added) to Spin Column CB3,

and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 s, then discard the flow-through and place the spin column into the collection tube.

2.1.7 Add 600 μ l Buffer PW (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 s. Discard the flow-through and place the spin column into the collection tube.

2.1.8 Repeat Step 7.

2.1.9 Centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 2 min to dry the membrane completely. Note: The residual ethanol of buffer PW may have some affection in downstream application

2.1.10 Place the Spin Column CB3 in a new clean 1.5 ml microcentrifuge tube, and pipet 50-200 μ l Buffer TE directly to the center of the membrane. Incubate at room temperature (15-25°C) for 2-5 min, and then centrifuge for 2 min at 12,000 rpm ($\sim 13,400 \times g$).

2.2 Bisulfite treatment of DNA

2.2.1 Add 130 μ l of the CT Conversion Reagent to 20 μ l of your DNA sample in a PCR tube. If the volume of the DNA sample is less than 20 μ l, make up the difference with water. Mix the sample by flicking the tube or pipetting the sample up and down, then centrifuge the liquid to the bottom of the tube.

2.2.2 Place the sample tube in a thermal cycler and perform the following steps*:

1. 98°C for 10 minutes
2. 64°C for 2.5 hours
3. 4°C storage up to 20 hours.

2.2.3 Add 600 μ l of M-Binding Buffer to a Zymo-Spin™ IC Column and place the column into a provided Collection Tube.

2.2.4 Load the sample (from Step 2) into the Zymo-Spin™ IC Column containing the M-Binding Buffer. Close the cap and mix by inverting the column several times.

2.2.5 Centrifuge at full speed ($>10,000 \times g$) for 30 seconds. Discard the flow-through.

2.2.6 Add 100 μ l of M-Wash Buffer to the column. Centrifuge at full speed for 30 seconds.

2.2.7 Add 200 μ l of M-Desulphonation Buffer to the column and let stand at room temperature (20-30°C) for 15-20 minutes. After the incubation, centrifuge at full speed for 30 seconds.

2.2.8 Add 200 μ l of M-Wash Buffer to the column. Centrifuge at full speed for 30 seconds. Add another 200 μ l of M-Wash Buffer and centrifuge for an additional 30 seconds.

2.2.9 Place the column into a 1.5 ml microcentrifuge tube. Add 10 μ l of M-Elution Buffer directly to the column matrix. Centrifuge for 30 seconds at full speed to elute the DNA.

2.3 PCR amplification

2.3.1 0.2ml PCR tube, 20ul pcr amplification system.

10×PCR Buffer	5.0 μ l
2.5mM dNTPs	5.0 μ l
10 μ M primers(F+R)	2.0 μ l+ 2.0 μ l

DNA	1.5 μ l
5U/ μ l Taq HS	0.2 μ l
ddH ₂ O	34.3 μ l

2.3.2 PCR amplification

Pre-denaturation	95°C, 5min	} 40×cycle
Denatured	95°C, 30s	
Annealing	55°C, 30s	
Extended	72°C, 30s	
End extension	72°C, 5min	

2.4 Electrophoresis: 1×TAE, 2.0% agarose, 4V/cm