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**Paraffin -fluorescence probe-FISH and Immunofluorescence protocol****1. Apparatus and reagents****1.1 Apparatus**

<b>Name</b>	<b>Producer</b>	<b>Model</b>
Dehydrator	WHJJ	JJ-12J
Paraffin embedding machine	WHJJ	JB-P5
Pathologic microtome	Leica	RM2016
Frozen flat	WHJJ	JB-L5
Water Bath-Slide	Kedee	KD-P
RNase-free glass microscope slides	Wanwu	
Micro-centrifuge	Wanwu	D1008E
Rocker	Wanwu	TSY-B
Vortex	Dragon	MX-F
Pipettor	Nikon	KE0003087/KA0056573
Microscopy Imaging system	Nikon	NIKON ECLIPSE CI NIKON DS-U3
Liquid blocker pen	Gene tech	GT1001
Refrigerator	HAIER	BCD-192TGN
Incubator	LABOTERY	GSP-70
Autoclave	PANASONIC	MLS-3751L-PC

**1.2 Major reagents**

<b>reagent</b>	<b>manufacturer</b>	<b>article number</b>
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
PBS solution	Wanwu	G0002
Proteinase K	Wanwu	G1205
4% of paraformaldehyde (DEPC water)	Wanwu	G1113
DEPC	Amresco	E174
20×SSC solution	Wanwu	G3016-4
BSA	Wanwu	G5001

---

DAPI	Wanwu	G1012
Anti-fluorescence sealing tablets	quenching Wanwu	G1401
hybridization buffer	Wanwu	G3016-3
First antibody		
Second antibody		

---

## 2. The steps of the experiment

2.1 **Organization fixation:** take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) above 12h.

2.2 **Dehydration:** the tissue is dehydrated by gradient alcohol, paraffin, embedding.

2.3 **Section:** the paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4 **Dewaxing and dehydration:** soak sections in 2 changes of xylene, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85%and 75% ethanol 5 minutes each. Wash in DEPC dilution.

2.5 **Digestion:** according to the tissue fixation time, the slices are boiled in the retrieval solution for 10-15 minutes and naturally cooled. Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, Add proteinase K(20ug/ml) working solution to cover objectives and incubate at 37°C for \_\_\_min. Wash in pure water, then wash three times in PBS (pH 7.4) on a Rocker device, 5 min each.

2.6 **Pre-hybridization:** add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

2.7 **Hybridization:** remove the pre-hybridization solution, add the probe hybridization solution with concentration of \_\_\_, and incubate the section in a humidity chamber and hybridize overnight at \_\_\_°C.

2.8 **Washing:** remove the hybridization solution. Wash sections in 2×SSC for 10 min at 37°C , wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.9 **Blocking:** add blocking serum to the section and incubate at room temperature for 30 min.

2.10 **Incubate first antibody:** PBS solution containing a \_\_\_dilution of primary antibody were added and incubated at 4°C overnight. Samples were then washed with PBS three times for 5 min each at RT.

2.11 **Incubate second antibody:** after washing, the section was incubated for 50 min with second antibody at RT. Samples were then washed with PBS three times for 5 min each at RT.

2.12 **Stain cell nuclei (counter stain):** incubate with DAPI for 8min in the dark, and then mounting.

2.13 **Microscopic examination and photography:** to take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

### 3. Interpretation of the results

The nuclear stained by DAPI were blue under ultraviolet excitation, and the positive expression was a kind of fluorescence labeled by corresponding luciferin. FAM (488) is green light, cy3 is red light. The results of mRNA in situ hybridization were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression , Different fluorescence brightness is strong or weak

**Attached table 1 probe information.**