Co-Immunoprecipitation 实验报告

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Co- Immunoprecipitation protocol**1、 Experiment instruments**

Name	Company	Model number
Microplate Reader	Rayto	RT-6100
Homogenizer	Wanwu	KZ-II
Centrifuge	DRAGONLAB	D3024R
Mini Centrifuge	Wanwu	D1008E
Vortex mixer	Wanwu	MX-F
Magnetic stirrers	Wanwu	MS-PB
Decolorizing shaker	Wanwu	TSY-B
Vertical Electrophoresis Apparatus	Wanwu	BV-2
Transfer Electrophoresis Apparatus	Wanwu	BT-2
Camera obscura	Wanwu	WGA0017
Darkroom lamp	Wanwu	WGA0016
Photographic film	Wanwu	WGJP0001
Ice Maker	SIMAG	SPR80
Sealer	Wanwu	WGXYQ0004
Scanister	EPSON	V370
Image-analysis software	Adobe	Adobe PhotoShop
Gray analysis software	Alpha Innotech	alphaEaseFC

2、 Reagents

Reagents	Company	Catlog
IP Lysis Buffer	Wanwu	G2038-100ML
Protein A/G-beads	Millipore	IP05
50*cocktail	Wanwu	G2006
PMSF (100mM)	Wanwu	G2008
Phosphatase Inhibitor	Wanwu	G2007
BCA protein quantitative detection kit	Wanwu	G2026
5*/2*loading buffer	Wanwu	G2013/G2031
SDS-PAGE kits	Wanwu	G2003
Protein Marker	Wanwu	26617

PVDF membrane 0.45μm	Wanwu	G6015-0.45
Milk	Wanwu	G5002-100G
TWEEN 20	Wanwu	WGT8220
ECL	Wanwu	G2014
Developer and fixer	Wanwu	G2019
PeroXidase-conjugated Goat Anti-Rabbit IgG(H+L)	Wanwu	GB23303
PeroXidase-conjugated Donkey Anti-Goat IgG(H+L)	Wanwu	GB23404
PeroXidase-conjugated Goat Anti-Mouse IgG(H+L)	Wanwu	GB23301
PeroXidase-conjugated Goat Anti-Rat IgG(H+L)	Wanwu	GB23302
Transfer buffer	Wanwu	G2017
Electrophoretic Buffer	Wanwu	G2018
TBS buffer solution	Wanwu	G0001

3、Protocol For CoIP

(1) Cell protein extraction

1.1 Cell culture (Adherent cells: The cells should be about 80% - 90% confluent in a 100 mm cell culture dish, or the number of cells in the culture bottle was $2-5 \times 10^7$. Suspension cells: Centrifuge to collect cells ($2-5 \times 10^7$)).

1.2 Add enough cold IP buffer (Protease and phosphatase inhibitors should be added before use). Maintain constant agitation for 10 min in ice bath, use pipette to mix it up. Then gently transfer the cell suspension into a 1.5ml microcentrifuge tube. Continue splitting for 20 min in ice bath.

1.3 Centrifuge at 12,000 rpm, 4°C for 10 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice; then determination with the BCA kit

1.4 Denature some lysis solution for Input

(2) Tissue protein extraction

2.1 Wash the sample with cold PBS for twice or three times, then homogenize with an electric homogenizer. Volumes of lysis buffer must be determined in relation to the amount of tissue present.

2.2 Gently transfer the cell suspension into a 1.5ml microcentrifuge tube. Maintain constant agitation for 30 min in ice bath, use pipette to mix it up.

2.3 Centrifuge at 12,000 rpm, 4°C for 10 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice; Use BCA method to determine protein concentration.

2.4 Denature some lysis solution for Input.

(3) Protocol For CoIP

- 3.1 Add 1.0 μ g IgG and 20 μ L Protein A/G PLUS-agarose to the IgG group. Add 20 μ L Protein A/G PLUS-agarose to the experimental group. Shaken at 4 °C for 1 h..
- 3.2 Centrifuge at 2000g ,4°C for 5 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant .
- 3.3 Incubate 1-10 μ L (0.2-2 μ g) antibody coupled beads in 250 μ L cell lysis supernatants, rotate overnight at 4 °C
- 3.4 Mix the lysates with 80 μ L Protein A/G PLUS-agarose rotating at 4 °C for 2 h.
- 3.5 Centrifuge at 2000g, 4°C for 5 min. gently remove the supernatant.
- 3.6 Wash the immunoprecipitation complex with 1 ml of precooled IP lysate (without any inhibitors) for 4 times. Centrifuge at 2000g ,4°C for 10 min , and the supernatant should be carefully discarded each time.
- 3.7 After the last washing ,discard the supernatant as much as possible, then add 80 μ L 1 $\times$ reduced loading buffer, boil for 10 min in boiling water. Centrifuge at 1000 g, 4 °C, for 5 min. then aspirate the supernatant.

4、Detection

- 4.1 Separate the protein sample by SDS-PAGE (5-15%) and analyze by western blot.
- 4.2 The results were analyzed in the form of Excel.