Immunoprecipitation 实验报告**1、实验器材**

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Immunoprecipitation protocol

1、Experiment instruments

Name	Company	Model number
Microplate Reader	Rayto	RT-6100
Homogenizer	Wanwu	KZ-II
Centrifuge	DRAGONLAB	D3024R
Mini Centrifuge	Wanwu	D1008E
Vortex mixer	Wanwu	MX-F
Magnetic stirrers	Wanwu	MS-PB
Table concentrator	Wanwu	TSY-B
Vertical Electrophoresis Apparatus	Wanwu	BV-2
Transfer Electrophoresis Apparatus	Wanwu	BT-2
Camera obscura	Wanwu	WGA0017
Darkroom lamp	Wanwu	WGA0016
Photographic film	Wanwu	WGJP0001
Ice Maker	SIMAG	SPR80
Sealer	Wanwu	WGXYQ0004
Image-analysis software	Adobe	Adobe PhotoShop
Gray analysis software	Alpha Innotech	alphaEaseFC

2、Reagents

Reagents	Company	Catlog
IP Lysis Buffer	Wanwu	G2038-100ML
Protein A/G-beads	Millipore	IP05
50*cocktail	Wanwu	G2006
PMSF (100mM)	Wanwu	G2008
Phosphatase Inhibitor	Wanwu	G2007
BCA protein quantitative detection kit	Wanwu	G2026
5*/2*loading buffer	Wanwu	G2013/G2031
SDS-PAGE kits	Wanwu	G2003
Protein Marker	Wanwu	26617
PVDF membrane 0.45μm	Wanwu	G6015-0.45
Milk	Wanwu	G5002-100G

TWEEN 20	Wanwu	WGT8220
ECL	Wanwu	G2014
Developer and fixer	Wanwu	G2019
PeroXidase-conjugated Goat Anti-Rabbit IgG(H+L)	Wanwu	GB23303
PeroXidase-conjugated Donkey Anti-Goat IgG(H+L)	Wanwu	GB23404
PeroXidase-conjugated Goat Anti-Mouse IgG(H+L)	Wanwu	GB23301
PeroXidase-conjugated Goat Anti-Rat IgG(H+L)	Wanwu	GB23302
Transfer buffer	Wanwu	G2017
Electrophoretic Buffer	Wanwu	G2018
TBS buffer solution	Wanwu	G0001

3、IP Protocol

(1) Preparation of lysate from cells

1.1 Cell culture (Adherent cells: The cells should be about 80% - 90% confluent in a 100 mm cell culture dish, or the number of cells in the culture bottle was $2-5 \times 10^7$, Suspension cells : Centrifuge to collect cells ($2-5 \times 10^7$))

1.2 Add enough cold IP buffer. (Protease and phosphatase inhibitors should be added before use) , vibrate. Maintain constant agitation for 10 min in ice bath, use pipette to mix it up. Then gently transfer the cell suspension into a 1.5ml microcentrifuge tube. Continue splitting for 20 min in ice bath.

1.3 Centrifuge at 12,000 rpm, 4°C for 10 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice; then determination with the BCA kit

1.4 Denature some lysis solution for Input

(2) Tissue protein extraction

2.1 Wash the sample with cold PBS for twice or three times , then homogenize with an electric homogenizer. Volumes of lysis buffer must be determined in relation to the amount of tissue present.

2.2 Gently transfer the cell suspension into a 1.5ml microcentrifuge tube. Maintain constant agitation for 30 min in ice bath, use pipette to mix it up.

2.3 Centrifuge at 12,000 rpm ,4°C for 10 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice; Use BCA method to determine protein concentration.

2.4 Denature some lysis solution for Input.

(3) Protocol For IP

- 3.1 Add 1.0 μ g IgG and 20 μ L Protein A/G PLUS-agarose to the IgG group. Add 20 μ L Protein A/G PLUS-agarose to the experimental group. Shaken at 4 °C for 1 h.
- 3.2 Centrifuge at 2000g ,4°C for 5 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant
- 3.3 Incubate 1-10 μ L (0.2-2 μ g) antibody coupled beads in 250 μ L cell lysis supernatants, rotate overnight at 4 °C
- 3.4 Mix the lysates with 80 μ L Protein A/G PLUS-agarose rotating at 4 °C for 2 h.
- 3.5 Centrifuge at 2000g, 4°C for 5 min. gently remove the supernatant.
- 3.6 Wash the immunoprecipitation complex with 1 ml of precooled IP lysate (without any inhibitors) for 4 times. Centrifuge at 2000g ,4°C for 10 min , and the supernatant should be carefully discarded each time
- 3.7 After the last washing ,try best to remove the supernatant. Elute protein with 80 μ L 1x SDS protein gel sample loading buffer, boil for 10 min in boiling water. Centrifuge at 1000 g, 4 °C, for 5 min. then aspirate the supernatant.

4、Detection

- 4.1 Load 50 μ L protein sample into the wells of the SDS-PAGE gel, along with molecular weight marker.
- 4.2 Determine those sample by the dying method with coomassie brilliant blue. Extraction of protein from gel for mass spectrometry